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40 minutes RT-qPCR Assay for Screening Spike N501Y and HV69-70del Mutations

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(1) WHO National Influenza Center and Virology Reference Laboratory of Ministry of Health Turkey, Ankara, Turkey

(2) Bioeksen R&D Technologies Limited, Istanbul Technical University Teknokent, Istanbul, Turkey

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Table 4. Instruments and their software validated with the Bio-Speedy® SARS-CoV-2 N501Y Mutation Detection Kit

Manufacturer	Instrument	Instrument Catalog No	Software	Software Catalog No
Roche	LightCycler® 96	5815916001	LightCycler® 96 Application and Instrument Software	Included in 5815916001
Bio-Rad	CFX96 Touch™	1855195	CFX Maestro™ Software	12004110
Qiagen	Rotor-Gene® Q 5 Plex Platform	9001570	Q-Rex Software v1.0	Included in 9001570
Bio Molecular System	Magnetic Induction Cyclers (Mic)	MIC-4	micPCR Software Version 2.8	Included in MIC-4
HiMedia	Insta Q96® Plus	LA1073	PCR Insta Q96 v1.0	Included in LA1073

Table 5. Reaction set-up and the details of the qPCR program

Reaction Setup			qPCR Program		
Component	Reaction N501	Reaction Y501	Cycle	Temperature	Duration
2X Prime Script Mix	5 µL	5 µL	1	52 °C	5 min
N501 Oligo Mix	2.5 µL				
Y501 Oligo Mix		2.5 µL	1	95 °C	10 sec
				95 °C	1 sec
Template Nucleic Acid	2.5 µL	2.5 µL	40	60 °C	1 sec
Total reaction volume	10 µL	10 µL		FAM/ROX/CY5 Read	