


RNA 保存液「RNA SHIELDER」



The graphic features a light blue background with a white speech bubble containing the text "いいえ違います!" (No, it's different!). To the right, a pink question mark icon is positioned above the text "組織中RNAを保存するもの?" (Something to preserve RNA in tissue?). The main title "RNA SHIELDER" is prominently displayed in yellow and black, with "RNA保存液" (RNA preservation liquid) written in smaller black text above it. Below the title, the text "RNA溶液を長期保存するための保存液です。" (It is a preservation liquid for long-term storage of RNA solutions.) is underlined in yellow.

RNA SHIELDER は、
Total RNA や ssRNA などの RNA 溶液を、冷蔵(2~8℃)で安定的に長期保存することができます。

- ▶凍結融解はもう必要ありません！
精度管理用コントロール RNA も冷蔵保存できます。
- ▶長期保存が実現！
従来は -20℃凍結保存でも分解され、長期保存は難しいとされてきました。
- ▶輸送にも便利！
冷蔵輸送(2~8℃)等によるダメージも最小限に抑えます。
- ▶逆転写反応や PCR を阻害しません。
逆転写 PCR やリアルタイム PCR、次世代シーケンス実験系等でも、そのまま使用することができます。

詳しくは[こちら](#)